


UNIVERSITY OF DELAWARE MASS SPECTROMETRY FACILITY

On-filter in-cell (OFIC) digestion of *intact C. elegans* with E4technology

Summary

This SOP pertains to processing biological samples for high resolution mass spectrometry (MS) based proteomics analysis. The protocol describes digesting proteins directly in the intact *C. elegans*, without cell lysis and protein extraction, thus it significantly simplifies the entire workflow, and drastically reduces sample loss. This method could be used for low-input, and single-worm/single-cell proteomics.

Note: ready-to-go E4 filter devices are available [here](#).



1. Worm treatment

Rinse live *C. elegans* three times with M9 buffer, then one time with water, to remove *E. coli* food source.

2. Sample loading

Pick worm manually, and transfer directly into E4filters that are pre-filled with 200 μ l of pure methanol. Visually inspect worm pick to confirm transfer.

Note: Estimated capacity for E4tip, < 100 worms, and E4 spin columns, 100-2,000 worms.

3. Worm fixation

Incubate E4filters at room temperature for 15 min. Centrifuge at 1,500 x g for two min, discard flow through. Add 200 μ l of methanol, repeat this step one more time.

Note: Flow through may be collected here for metabolomics analysis.

4. Reduction and alkylation

Add 100 μ l of 50 mM triethylammonium bicarbonate (TEAB) and final concentration of 10 mM Tris(2-carboxyethyl)phosphine (TCEP) and 40mM chloroacetamide (CAA), incubate at 45°C for 10 min with gentle shaking.

5. Wash

Add 200 μ l of 50 mM TEAB solution, centrifuge at 1,500 x g for two min, discard flow through.

6. Digestion

Add 100 μ l 50 mM TEAB, desired enzyme (Trypsin or Trypsin/Lys-C mix) at 1:50 ratio. Incubate at 37°C for 16-18 hours with gentle shaking.

Note: no cap is required for E4tips.

7. Acidification and desalting

After digestion, add formic acid to final concentration of 1%, centrifuge at 500 x g for 10 min. Add 200 μ l 0.5% acetic acid in water, centrifuge at 1,500 x g for 2 min, discard flow through.

Note: here, E4tips (E4tip ultra) can be transferred to Evosep ONE LC for direct LCMS acquisition.

8. Elution

Transfer E4filters to clean collection tubes, do two sequential elution by adding 200 μ l 60% acetonitrile/0.5% acetic acid in water (elution I), and 80% acetonitrile/0.5% acetic acid in water (elution II), centrifuge at 1,500 x g for 2 min to collect elution to the same tube. Dry samples in SpeedVac, and store at -80°C. The peptides are now clean and ready for LCMS analysis.

