AVIII 400 NMR User Guide

November 11, 2012

The CBC NMR Laboratory University of Delaware

Contents

1	INTRODUCTION	.2
1.1 1.2 1.3	SAFETY NMR SAMPLE PREPARATION	.2 .2 .2
1.4 1.5		.3 .3
2	TOPSPIN NMR SOFTWARE	.4
2.1 2.2	TOPSPIN INTERFACE	.4 .5
3	STEP-BY-STEP GUIDE FOR ¹ H NMR	.9
4	STEP-BY-STEP GUIDE FOR ¹³ C SPECTRUM	17

List of Figures

Figure 1: TopSpin Interface
Figure 2: Lock Display
Figure 3: BSMS Panel
Figure 4: Fill in the "New" Window for First Time
Figure 5: User Preferences – Part I
Figure 6: User Preferences – Part II
Figure 7: Status Bar Preferences
Figure 8: Lock Display Preferences
Figure 9: BSMS Display Preferences
Figure 10: Status Bar
Figure 11: The "New dataset" Window
Figure 12: The "rpar" Window
Figure 13: Sample Insertion
Figure 14: Sample Icon in Status Bar10
Figure 15: Solvent Table
Figure 16: Lock Display Before and After Lock
Figure 17: Probe Tuning Curve
Figure 18: Data Acquisition in Progress
Figure 19: Processed 1H NMR spectrum
Figure 20: Click and Drag a Spectral Region
Figure 21: Click and Drag to Define a Region of Peak Picking
Figure 22: Spectrum with Peak Pickings
Figure 23: Spectrum Integration
Figure 24: Calibration of Integral
Figure 25: Processed ¹ H Spectrum
Figure 26: Printing Spectrum
Figure 27: Exit TopSpin
Figure 28: Create a New Dataset for ¹³ C Spectrum
Figure 29: Probe Tuning and Matching for ¹³ C and ¹ H Frequencies
Figure 30: Sample Icon in Status Bar
Figure 31: Acquisition Window When Spectrum Button is Highlighted
Figure 32: Processed ¹³ C Spectrum

PREFACE

Purpose of This Document

This user guide aims to familiarize you with data acquisition using TopSpin software on the AVIII 400 NMR Spectrometer. The document does not cover all the basic information and details of the application.

Intended Audience

This document is intended for users of the AVIII 400 NMR spectrometer who are familiar with the basic principles of the spectrometer operation. It will enable you to understand the details of basic operations. This document is *not* intended to replace any form of official Bruker user guides or related Bruker documents.

1 Introduction

The AVIII400 NMR spectrometer is equipped with a Bruker BBFO probe. This instrument is used for analyzing time and air sensitive samples. It is also suitable for the study of reaction kinetics of chemical reactions and for a variable temperature (VT) NMR experiment. Prior user training is required to operate the spectrometer. Additional training is available for the VT experiment. The web reservation is required to use the AVIII 400 NMR spectrometer. The spectrometer is running on a Bruker TopSpin software.

1.1 Safety

Iron and other ferro magnetic objects must **NOT** be brought into the vicinity of the magnets. The strong magnetic fields may erase the information of credit card, student ID card, and other magnetic media. No admission is allowed for persons with pacemakers and other metallic implants.

1.2 NMR Sample Preparation

- Always use clean and dry sample tubes
- Use medium to high quality sample tubes
- Filter sample solution if particles are present
- Keep the sample volume approximately 0.5 0.6 mL.
- Wipe the sample tube clean before inserting into magnet
- The sample tube should sit tightly inside the spinner
- Use a sample depth gauge to adjust position of sample spinner

1.3 NMR Probe

The AVIII 400 MHz NMR spectrometer is equipped with a BBFO probe. The BBFO probe has two radio frequency channels (¹H and X-nucleus) with built-in capability of automatic tuning and match (ATM). The high frequency channel is tuned for ¹H observe and decoupling. The observe channel (X) covers a frequency range from ¹⁰⁹Ag to ¹⁹F. The probe temperature range is from -150 to 150 °C. The commonly used acquisition files for this probe are prefixed with "*bbfo*"

1.4 Login

Once NMR training has been completed, username and password will be provided by the NMR spectroscopist. In general, the username is the same as the university email address, therefore, one must first obtain a university email address before an NMR account is created.

1.5 Workstation

The NMR workstation OS is a Linux system. Unlike Windows, please do not use special characters such as *, space, \$, %, etc. in the filenames. To separate concepts, one may consider using a dash or an underscore.

2 TopSpin NMR Software

2.1 TopSpin Interface

To start the TopSpin (TS) NMR software, click on the TopSpin icon on the desktop. A TopSpin interface (Figure 1), the lock display window (Figure 2), and the BSMS panel (Figure 3) will appear.

<u>e</u>	Bruker TopSpin 3.1.	b.58 on aviil400.chem.	udel.edu as sbai			
<u>S</u> tart <u>A</u> cquire <u>P</u> roce	ess A <u>n</u> alyse P <u>u</u> blis	sh <u>V</u> iew <u>M</u> ar	nage 🕜	Menu Bar		T 🗖 B
C <u>r</u> eate Data	set [📓 Find Dataset 🔄	Open <u>D</u> ataset 📘 P	aste Dataset	🖹 Read Pars.	Suhmanu Traa	
			∼≞ ₹₹∕∕	TS EFP	Spectral Tools	
Browser Last50 Groups 1 first 2	1 /home/sbai/mynmr					r 🛛 🖂
← 🗖 ColonialMetals	n ProcPars AcquPars Title	PulseProg Peaks	Integrals Sample	Structure Plot	Fid Acqu	
P □ /home/sbai/chem603 demo						
← 📑 /home/sbai/demo						
e 🖬 mseg		Spectrum Panel				
← 🗂 /home/sbai/nmrdata		Spectrum 1 and				
← ☐ /home/sbai/proj3 ← ☐ /home/sbai/project1						
←						
← □ proj3 ← □ projx						
← □ prox						
← □ /home/sbai/taber ← □ /home/sbai/traing		No r No proc	aw data available essed data availab	le		
Data Browsers						
Command Line						
Message Board						
Acquisition information	Fid Flash Lock Sample	Shim Coil POWCHK	Sample Tempe	rature Spooler	BSMS status messag	e Time 13:41:47
no acquisition running		302 K	On Beg	delayed:	Autoshim 📮 Locked 🚍	Nov 08

Figure 1: TopSpin Interface

The TopSpin interface contains several portions: the Menu Bar (on the top) with the buttons of **Start, Acquire, Process** etc. The Submenu Tree (the second row) shows special functions such as **Create Dataset, Find Dataset**, etc., corresponding to the **Start**, is located below the Menu Bar. Below the submenu, there are two rows of commonly used data manipulating tools. The left panel provides a number of ways to browse for data. The right panel (Spectrum panel) contains **Spectrum** tab, **ProcPars** (Processing parameter) tab, **AcquPars** (acquisition parameter) tab etc. Below the Spectrum panel, a command line is located. A command line allows user to enter a command manually. Below the command line, a message board provides real-time experiment information. The bottom section of the interface contains a number of status indicators, such as *Lock, Sample, Sample Temperature* etc.

	B BSMS Control Suite - 🗆 🗙
Lock Display	Main Lock/Level Shim Autoshim Service Log Help
▋▋▋▋▋	AUTO
	Lock Phase Power Gain Shim
	LOCK
	On-Off Phase Power Gain
	SAMPLE
	LIFT SPIN Measure Rate Lock Lost
	SHIM
	Spin. Z Z ² Z ³ Z ⁴
	X XZ
	Y YZ
	XY
	χ2-γ2
	STD BY
	Previous Actual Step
	Absolute + Reset
	Difference –
	Stepsize
	STD BY
	Config
	🗹 External
The sele second of the particular stands and a second from the second state and the	
	Sample: down missing up Shim coil temperature
	302 K

Figure 2: Lock Display

Figure 3: BSMS Panel

The lock display window shows ²H signal profile of deuterium frequency lock and the BSMS panel allows user to execute the sample related commands, such as sample insertion and ejection, manual shimming, locking, lock gain adjusting, etc.

2.2 Recommended Preferences

When TopSpin is started for the first time, the following preferences are recommended:

 In the command line, type *edc* and a window called "New..." opens. Fill in the Name field with a text such as "first". Assign both EXPNO and PROCNO to a numeric number such as "1". In the Experiment field, select a proton experiment acquisition file, such as "*bbfo1hstd*". The data directory field (DIR) specifies data storage location (/home/username/data_directory). As an example shown in Figure 4, the username is *sbai* and data_directory is *proj1*. *Note*: the directory path /home/username is fixed sequence as user's home directory, which cannot be changed. The data_directory, however, can be set by user. When all fields are filled in, click on OK to proceed.

۵.	Nev	V X
Prepare for a new experinitializing its NMR para For multi-receiver exper Please define the numb	iment by creating meters according t riments several dat er <u>of rece</u> ivers in th	a new data set and to the selected experiment type. tasets are created. te Options.
	first	
EXPNO	1	
PROCNO	1	
TITLE		
O Use current paramete	ers	
Experiment bbfo1hs	td	Select
Options		
🗹 Set solvent:		CDCI3
○ Execute "getpro	sol"	
○ Keep parameter	s:	P 1, O1, PLW 1 💌 Change
DIR		/home/sbai/proj1
🗷 Show new datas	et in new window	
Receivers (1,2,1	.6)	1
	<u></u> K	<u>C</u> ancel More <u>I</u> nfo <u>H</u> elp

Figure 4: Fill in the "New..." Window for First Time

2. Type *set* in the command line, the User Preference window opens. Make sure "*Auto-open last used dataset when restarting TopSpin*" is checked as shown in Figure 5

ώ.	User preferences	
Administration items Window settings Text editors Miscellaneous Remote connection Directories	Administration items Auto-open last used dataset when restarting TopSpin Show TopSpin data examples directory in data browser Setup users for TopSpin-internal login/logoff and esign Automatic termination of TopSpin when idle time exceeded	Change Change
Acquisition	Automatic locking of TopSpin when idle time exceeded	Change

Figure 5: User Preferences – Part I

3. Also in the User Preference window, under Acquisition, make sure four items (*"Show "ased" ...; "Overwrite..."*; *"Display Digital..."*; and *"Auto Open Acquisition..."* are checked as shown in Figure 6.



Figure 6: User Preferences – Part II

4. Then click on **Change** associated with the **Status Bar Preferences**, a window opens as shown in Figure 7. Make sure the items in Figure 7 are checked. Click on **Apply** then on **Back**.

🖕 Status Bar Preferences	X
Acquisition status bar	
Auto open acquisition status bar	
Include spooler	
Include time	\mathbf{r}
Include sample temperature	\mathbf{r}
Include channel or accessory chanel temperature	
Include acquisition status	~
Include acquisition indicator	~
Include lock signal	~
Include MAS spinning rate	
Include peak power check (POWCHK) indicator	
Include sample state	~
Include shim coil temperature	~
Include BSMS status	
Include amplifier control	
Apply Back Clos	e

Figure 7: Status Bar Preferences

5. Then click on **Change** associated with **Lock Display Preferences**. In the pop-up window (Figure 8), make sure "*Auto open LOCK display*" is checked. Click on **Apply** and then on **Back**.

🍯 🛛 Lock Display Preferen	ces 🛛 🗙
Lock display	
Auto open LOCK display	
Background color	Change
Color 1	Change
Color 2	Change
Color of grid lines	Change
Number of horizontal grid li	nes 4
Number of vertical grid lines	5
<u>A</u> pply <u>B</u> ack	<u>C</u> lose

Figure 8: Lock Display Preferences

6. Finally click on **Change** button of BSMS Display Preferences, check "auto open BSMS display" and "External BSMS display". Click on **Apply** and then on **Back**.

BSMS Display Preferences	×
<mark>BSMS display</mark> Auto open BSMS display External BSMS display (newer BSMS type	es only) 🗹
<u>A</u> pply <u>B</u> ack	<u>C</u> lose

Figure 9: BSMS Display Preferences

- 7. Back to the User Preference window, click on **Apply** (*c.f.* Figure 6) and close the window.
- 8. Pointing the cursor to the bottom portion of TopSpin Interface, right click the mouse, and select "*Show Status Bar*". The status bar, including *Acquisition Information, Lock, Sample, Sample temperature* etc. as shown in Figure 10, now open.

Acquisition information	Fid Flash	Lock	Sample	Shim Coil	POWCHK	Sample Temperature	Spooler	BSMS status message	Time
no acquisition running			1 ?	Temperature 302 K	\checkmark	Corr. 298.4 K	queued: 0 delayed: 0 cron: 0	Autoshim 😑 Locked 🝣 Error	16:46:33 Nov 13

Figure 10: Status Bar

3 Step-by-step Guide for ¹H NMR

Follow these steps to obtain a ¹H NMR spectrum:

1. Click on **Start** in the menu bar. Then click on **Create Dataset** button in its submenu bar. A "**New...**" window opens up.

Č	New	×
Prepare for a new experime initializing its NMR paramet For multi-receiver experime Please define the number of	nt by creating a new data set and ers according to the selected experiment type, ents several datasets are created. Freceivers in the Options.	
NAME	demo	
EXPNO	1	
PROCNO	1	
TITLE	demo	
○ Use current parameters		
Experiment bbfo1hstd	Select	
 Options 		
🗹 Set solvent:	CDCI3 🗸	
○ Execute "getprosol"		
⊖ Keep parameters:	P 1, O1, PLW 1 💌 Change	
DIR	/home/sbai/mynmr	•
🗹 Show new dataset ir	n new window	_
Receivers (1,2,16)	1	
	<u>Q</u> K <u>C</u> ancel More <u>I</u> nfo <u>H</u> e	lp

Figure 11: The "New dataset" Window

- 2. Fill in the fields of NAME, EXPNO, PROCNO, and TITLE. *Note: that only numerical numbers are allowed in the fields of EXPNO and PROCNO.*
- 3. Click on Select associated with the Experiment field, a "rpar" window opens:

4		Parameter Sets:	rpar					
<u>F</u> ile <u>O</u> ptions <u>H</u>	elp	Sou	rce = /opt/topspin3.1.b.	58/exp/stan/nmr/par/user				
Find file names 👻	enter any string, *, ?	Clear						
Class = Any - Din	Class = Any v Dim = Any v Show Recommended							
Type = Any V SubType = Any V SubTypeB = Any V Reset Filters								
bbfo13cstdIm	bbfo13cstdsm	bbfo19f	bbfo19fcosy	bbfo19fcosygp				
bbfo19fdec	bbfo1dnoe	bbfo1hpara	bbfo1hstd	bbfo29sistd				
bbfo2hlsstd	bbfo2Hobs_BB	bbfo31p	bbfo31pdec	bbfo77sestd				
bbfohsgcstd	bbfoto csystd	f19-cosyph	gracec13	gradtest				
noesygppr1d	user.tar.gz_FILES							

Figure 12: The "rpar..." Window

- 4. Highlight the acquisition file "*bbfo1hstd*" and then click on **Set Selected Item in Editor**.
- 5. Go back to "**New...**" window, select appropriate solvent from the '*solvent*' pulldown menu, and make sure '**getprosol**' is checked. Click on **OK** button to finish dataset creation.
- 6. Click on Acquire in the menu bar, select Turn on Sample Lift Air by clicking on the green triangle sign of the Sample button in the submenu tree. This will turn on the lifting air for sample insertion. Alternatively, type *ej* in the command line to turn on the lifting air. When lifting air is on, place the NMR sample with spinner on the top of magnet. Click on Sample in the submenu again and select Turn off Sample Lift Air. Or simply type "*ij*" in the command line. The sample will be lowered into the probe. *Note: in the bottom section of the TS interface, Sample icon is changed as indicated in* Figure 14.



Figure 13: Sample Insertion



Figure 14: Sample Icon in Status Bar

7. Click on Lock in the Acquire submenu. A Solvent Table will appear. Select appropriate solvent, and the click on the "OK" button.

۵	Solvents table	×
▲ Solvent	Description	
Acetic	acetic acid-d4	
Acetone	acetone-d6	
C6D6	benzene-d6	
CD2CI2	Methylene Chloride d2	
CD3CN	acetonitrile-d3	
CD3CN_SPE	LC-SPE Solvent (Acetonitrile)	
CD3OD_SPE	LC-SPE Solvent (Methanol-d4)	
CDCI3	chloroform-d	
CH3CN+D2O	HPLC Solvent (Acetonitril/D2O)	
CH3OH+D2O	HPLC Solvent (Methanol/D2O)	
D20	deuteriumoxide	
DMF	dimethylformamide-d7	
DMSO	dimethylsulfoxide-d6	
Ft OD	athanal_d6	

Figure 15: Solvent Table

8. Locking signal is then changed as indicated in Figure 16.



Figure 16: Lock Display Before and After Lock

9. Click on **Tune** in the **Acquire** submenu tree to perform a probe automatic tuning and matching. Or type *atma* in the command line. When the probe tuning starts, a tuning curve (Figure 17) appears and soon disappears in **Acqu** panel.



Figure 17: Probe Tuning Curve

- 10. After the probe tuning, click on **Shim** (or type *topshim* in the command line) to do an automatic shimming which may take 30 seconds or longer.
- 11. Click on **Gain** (or type *rga* in the command line) to automatically adjust receiver gain.

- 12. When receiver gain adjustment is finished, click on **Go** (or type *zg* in command line) to acquire proton NMR data. *Note: the status bar shows the progress of data acquisition*.
- 13. When data acquisition is in progress, the **Acqu** Tab of the spectrum panel is activated and shows the FID data or spectrum data depending on which button is depressed as shown in Figure 18.



Figure 18: Data Acquisition in Progress

14. After data acquisition is finished, click on **Spectrum** tab in the spectrum panel. A message of "*1D raw data available and no processed data available*" is shown.

Click on **Process** in the menu bar followed by clicking on **Proc. Spectrum** in **Process** submenu tree. A processed spectrum is shown in **Spectrum** tab.

Q	Bruker TopSpin 3.1.b,58 on aviii400.chem.udel.edu as sbai 🖉 🖬					
<u>S</u> tart <u>A</u> cquire	<u>P</u> rocess <u>An</u> alyse	P <u>u</u> blish <u>V</u> iew	<u>M</u> anage 🕜		1 2 L B	
 ∧ Pro <u>c</u> . Spectru	ım 🗢 🐴 Adjust Phase 🤜	🔻 🕅 Calib. A <u>x</u> is 👻 🎌 P	ck P <u>e</u> aks ▼ ∫ <u>I</u> ntegrate ▼	A <u>d</u> vanced ▼		
Image: Second	©кн © (+→\$) ⊋ _⊖ ∰ <mark>≞</mark> (+→\$)					
Browser Last50 Groups	2 Acquisition finished: /home/	sbai/mynmr/demo/1/pdata/1			r 0' 🛛	
Ŷ □ mynmr	Spectrum ProcPars Acq	uPars Title PulseProg Pea	ks Integrals Sample Structure	Plot Fid Acqu		
 background cCl demo first Fox se.test decond decond decond vc second vc vt 	demo 0,36 pap / 3747.17 Hz Index = 22861 - 22425 Value = -0.0001505 cm					
			يه المستبالية المالية المالية	- Cleman		
	15	10	5	0	[ppm]	
L						
Acquisition inform	ation Fid Flash Loci	Sample Shim Coil Temperature 303 K	MCHK Sample Temperature	Spooler BSMS status messa Jeued: 0 Llayed: 0 on: 0 Autoshim C Locked C	ge Time 14:01:08 Nov 08	

Figure 19: Processed 1H NMR spectrum

15. Use the left button of the mouse to click and drag a region that covers all resonances, and then release the button. The spectrum will be expanded as defined by vertical cursors.

Q	Bruker TopSpin 3.1.b.58 on aviii400.chem.udel.edu as shai					
<u><u>S</u>tart <u>A</u>cquire</u>	Process Analyse Publish View Manage 🕜	1 2 B				
A Pro <u>c</u> . Spectre	um 👻 🕎 Adjust Phase 👻 🔥 Calib. A <u>x</u> is 👻 🎊 Pick P <u>e</u> aks 👻 🥤 Integrate 👻 A <u>d</u> vanced 🛩					
Image: Second secon						
Browser Last50 Groups	2 Acquisition finished: /home/sbai/mynmr/demo/1/pdata/1	r ø'⊠				
የ ⊑ mynmr	Spectrum ProcPars AcquPars Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu					
Cackground	demo					
← 🗖 demo	-0.47 ppn / -188.54 Hz Index = 54550 - 54614					
C first C Fox	Value = 0.01056 cm					
← 🗖 se_test						
← 🖾 second ← 🗂 topshim						
∾ 🗂 vt						
	Distance = 8.46 ppm / 3385.81 Hz					
	15 10 5 0	rbburl				
L						
Acquisition inform	nation Fid Flash Lock Sample Shim Coil POWCHK Sample Temperature Spooler BSMS status message	Time				
no acquisition running		14:07:17 Nov 08				
	j j katoshim Gocked Ger	ror				

Figure 20: Click and Drag a Spectral Region

16. Click **Peak Picking** in **Process** submenu tree, and make sure the selection icon (circled in Figure 21) in the **Peak Picking** submenu is highlighted. Then click and drag a green box, which defines the minimum and the maximum limits of peak picking, as shown in Figure 21.

(Bruke	er TopSpin 3.1.b.58 on a	vIII400.chem.udel.edu as sb	al	×
<u>S</u> tart <u>A</u> cquire	<u>P</u> rocess A <u>n</u> aly	se P <u>u</u> blish <u>V</u>	iew <u>M</u> anage 🕜		1 Z 🗖 B
Λ Pro <u>c</u> . Spect	rum 🗢 🔥 Adjust Phas	e 🗢 🕅 🚴 Calib. A <u>x</u> is 🚽	Pick P <u>e</u> aks → ∫ I	ntegrate 🗢 🛛 A <u>d</u> vanced 🗢	
Image: Second secon	₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽ ₽			TS EFP	
Browser Last50 Groups	2 Acquisition finished: /hon	ne/sbai/mynmr/demo/1/pd	ata/1		ة 10 N
P □ mynmr	<mark>╷╖(⊢)</mark> р⋇∣╨ _ѿ т	⋭⋇⊠∣DĘ(Ę))			
 CCEI Gemo Grist Fox Ge test 	demo				
 ■ second ■ topshim ■ vt 	0.307 ppm / 122.819 Hz Defire new region: Drag with left mouse butto Number of peaks: 0	n			
		6	4	2	0 [ppm]
	481	-	•	-	* (ppm)
-	((
Acquisition info no acquisition running	rmation Fid Flash L	ock Sample Shim C Tempera 303	iture K	nperature Spooler BSM 00.0 K Reg. State: ♥ cron: 0 Autosh	IS status message ▲ Y3 0 im ⊖ Locked ♥ Error Time 14:08:07 Nov 08

Figure 21: Click and Drag to Define a Region of Peak Picking

17. Click on **Save and Return** as circled in Figure 21 to save all peak pickings. A spectrum with peak pickings is shown in Figure 22.

le l	Bruker TopSpin 3.1.b.58 on avill400.chem.udel.edu as sbai	_ 🗆 🗙
📑 <u>S</u> tart <u>A</u> cquire	e Process Analyse Publish View Manage 🕢	1 2 L B
A Pro <u>c</u> . Spec	trum 👻 🕎 Adjust Phase 👻 🗼 Calib. Axis 👻 🎊 Pick Peaks 💌 🥤 Integrate 👻 Advanced 💌	
Image: Second secon	◙®™ <mark>®</mark> ⊨→? I ⊨ ⊗‰ ≥ ∵⊤ ∕ < ≈ ™	
Browser Last50 Groups	🖞 🔁 Acquisition finished: /home/sbai/mynmr/demo/1/pdata/1	ت ت N
P-III mynmr	<mark>╨┍</mark> ┉╳╨ѿѿᢂ┛┨╝╹	
background background c	demo Residence 7.88 part 298.33 Hz before her region: 0 ray wheth left souse button Waber of feaks: 52	60/M-
		······································
	6 4 2	0 [ppm]
l		
Acquisition info	ormation ng Fid Flash Lock Sample Shim Coll Temperature 303 K 303 K 0 0 0 Reg. State: 0 ron: 0 Autoshim Colcole	essage Time 14:08:31 Nov 08

Figure 22: Spectrum with Peak Pickings

14

18. To integrate spectrum, click on **Integrate** in the submenu tree of **Process**. Make sure the selection tool is highlighted (as circled in Figure 23), then click and drag an integral range around a peak to be integrated. Continue on the rest of the resonances.

6		Bruke	r TopSpin 3.1.b.	58 on aviii40	0.chem.udel.	edu as sbai		_ - ×
<u>S</u> tart	Acquire	Process A <u>n</u> alys	e P <u>u</u> blish	<u> </u>	<u>M</u> anage	0		12 B
۸ Pr	o <u>c</u> . Spectrum v	Adjust Phase	e ▼ 🚴 Calib.	A <u>x</u> is 🗢 🕅	Pick P <u>e</u> aks	⇒ ∫ įnteg	grate 🗢 🛛 A <u>d</u> vanced 🗢	
Browser Last50	Groups 2 A	equisition finished: /hom	ie/sbai/mynmr/dei	mo/1/pdata/1				o و ي
P ☐ mynmr		┖┛┝┛╫╘	s # 🚅 🛱 🖛	₽₽ ⊿ ∟	¥ *2 /2 ≎	\$ *X = ±	╀\$∦⊴┨┨╽	
 CCB demo first fi	der Hous Som DEFI DEFI Retu	770 e Sensitivity: 1.0 = 0.0000 ME RECON MODE ne: Orag USINg Tert mo rrn: Left-Click highlig	use button ted icon					o (ppm)
Acqu no acquisit	uisition information	h Fid Flash L	ock Sample	Shim Coil Temperature 303 K		ample Tempera Orr. 300. Reg. S	ture Spooler BSMS Q K queued: 0 delayed: 0 cron: 0 Autoshir	S status message Time 14:08:59 Nov 08 Nov 08

Figure 23: Spectrum Integration

19. To calibrate integral, right click over the integral on the screen and select **Calibrate**, fill in appropriate value in the pop-up box as shown in Figure 24. Finally click on '**Save and Return'** to save the integrations and return to spectrum panel.



Figure 24: Calibration of Integral

20. The final processed ¹H spectrum is shown in Figure 25.



Figure 25: Processed ¹H Spectrum

21. Finally, to print out the spectrum, click on **Publish** in the menu bar, and then click on **Print** in the submenu, as indicated as the red circled tool buttons below.

4		Bruker TopSp	in 3.1.h.58 on aviil	100.chem.udel.edu	as sbai			
🧾 <u>S</u> tart <u>/</u>	<u>A</u> cquire <u>P</u> roces	s A <u>n</u> alyse	P <u>u</u> blish <u>V</u> iew	<u>M</u> anage	2		1 🗖 B	
	Copy Print V Plot Layout V E-Mail							
Image: Second				▶ ■ᡂ₩₽	S REFP			
Browser Last50 Gro	ups 🚦 1 demo 1	1 /home/sbai/mynmr					r 🛛 🛛	
🕶 🗂 /home/sbai/proj1	Spectrum	ProcPars AcquPars	Title PulseProg	Peaks Integrals	Sample Structure	Plot Fid Acqu		
 ← ☐ /home/sbai/mynmr ← ☐ background ← ☐ CCEI ← ☐ demo 	demo	(

Figure 26: Printing Spectrum

- 22. In the BSMS panel, click on a green button "Lock ON-OFF" to turn off the frequency lock. Then click on a green button "Lift" to lift the sample from the probe. Remove the sample from the magnet, click on the Lift again to turn off the lifting air.
- 23. Exit the TopSpin program by clicking on the "x" sign as circled in Figure 27.

۵	Bruker TopSpin 3,1,b,58 on avIII400.chem.udel.edu as sbai									
	<u>S</u> tart	<u>A</u> cquire	<u>P</u> rocess	A <u>n</u> alyse	P <u>u</u> blish	<u>∨</u> iew	<u>M</u> anage	2		TB
			<u>[]</u> <u>С</u> ор	y 🗳 P <u>r</u> int →	P <u>l</u> ot	Layout 🗢	<u> </u>	📑 <u>E</u> -Mail		
	*8 * /8 /		₩N (Q) → N		Hz ppm		► ● ● ● ● ● ● ● ●	<u>۶</u> ۸	TS EFP	

Figure 27: Exit TopSpin

4 Step-by-step Guide for ¹³C Spectrum

We assume a proton NMR spectrum has been acquired by following the steps described in section 3. In this case, the frequency locking and shimming have been done. We can now follow these steps to obtain a 13 C spectrum.

 First, create a new dataset for ¹³C spectrum by first clicking on Start in the menu bar, followed by clicking on Create Dataset in its submenu tree. Change EXPNO from 1 to 2, and select an appropriate ¹³C acquisition parameter file (*bbfo13cstdsm*) in Experiment pull-down menu. Select appropriate solvent in Solvent pull-down window. Make sure Execute "getprosol" is checked. Click on OK to continue.

	New	×
Prepare for a new experimer initializing its NMR paramete For multi-receiver experimer Please define the number of	nt by creating a new data set and rs according to the selected experiment type nts several datasets are created. receivers in the Options.	
NAME	demo	
EXPNO	2	
PROCNO	1	
TITLE	demo	
O Use current parameters		
Experiment bbfo13cstdsn	n	Select
Options		
🗹 Set solvent:	CDCI3	
Execute "getprosol"		
⊖ Keep parameters:	P 1, O1, PLW 1 🔻 Chang	e
DIR	/home/sbai/mynmr	-
🗹 Show new dataset in	new window	
Receivers (1,2,16)	1	
	<u>O</u> K <u>C</u> ancel More <u>I</u> nfo.	<u>H</u> elp

Figure 28: Create a New Dataset for ¹³C Spectrum

2. Once a new dataset is created, click on **Acquire** in the menu bar and then on **Tune** in its submenu tree. Automatic tuning and matching starts momentarily. The tuning curves are shown for ¹³C frequency tuning and then for ¹H frequency sequentially (Figure 29).



Figure 29: Probe Tuning and Matching for ¹³C and ¹H Frequencies

3. Click on **Spin** in the **Acquire** submenu tree to spin the sample. Wait until **Sample** icon in the status bar changes as indicated in Figure 30.



Figure 30: Sample Icon in Status Bar

4. Click on **Go** in **Acquire** submenu tree (or type *zg* in the command line) to start the data acquisition. Make sure the *Spectrum* icon (circled in) is highlighted in **Acqu** panel. When signal-to-noise ratio reaches to a satisfactory level, type *halt* in the command line to pause the data acquisition. *Note: if* **Stop** *button is depressed, the experiment will be terminated without saving the data*.

ē.	Bruker To	pSpin 3.1.b.58 on aviii40).chem.udel.edu as sbai		_ = ×
<u>S</u> tart <u>A</u> cquire	<u>P</u> rocess A <u>n</u> alyse	P <u>u</u> blish <u>V</u> iew	<u>M</u> anage 🕜		1 2 3 A B
💐 Sample 👻 🌞 I	Lock V T <u>u</u> ne → 👪 Sp <u>i</u> r	n 🔻 👫 Shim 👻 🔏 Pr	osol 🗢 🗽 Gain 🗢 🕨 🕨	Go 🔻 Options 🔻	
Image: Second secon			Ìœ ₩₽ \$ ▼₹∕\$	TS EFP	
Browser Last50 Groups	ACQ: demo 2 1 /home/sbai/	mynmr			r 🛛 🖂
P ⊑ mynmr	Spectrum ProcPars Acqui	Pars Title PulseProg F	Peaks Integrals Sample	Structure Plot Fid Acqu	
← 🖬 background ← 🖬 CCEI	╞┉╞┉╞╱┊᠘᠊	. 💿 🗉			
	demo				
I = 2g30 - demo I = 2ggg30 - demo	157.2 ppm / 15813.6 Hz Index = 273453 = 274481				Acquisition
← ☐ first	Value = 1.764 Cn				PULPROG = zgpg30
← ☐ Fox ← ☐ se_test					NUC1 = 13C
← 🗖 second					SWH = 24038.5
← 🖬 topsnim ← 🗂 vt					SCANS = 26/256
					RES. TINE = 13n5s/14n50s
	200	120	100	50	0 [ppm]
L					
Acquisition inform	nation Fid Flash Lock	Sample Shim Coil I	OWCHK Sample Temperat	ure Spooler BSMS s	itatus message Time
Mame/Expno: demo/2 Scan: 26/256		Temperature	/ Corr. 300.0	OK queued: 0 delayed: 0	Y3 0 14:25:47
Residual time: 13m2s			V On 🛇 Reg. St	ate: 🕑 cron: 0 Autoshim	Cocked C Error

Figure 31: Acquisition Window When Spectrum Button is Highlighted

5. After the experiment is terminated, click on **Process** in the menu bar, followed by clicking on **Proc. Spectrum** in its submenu tree. A processed ¹³C spectrum is shown in Figure 32.

Bruker TopSpin 3.1.5.58 on aviil400,chem.udel.edu as sbai							
🧃 <u>S</u> tart <u>A</u> cquire	Process A <u>n</u> alyse	P <u>u</u> blish <u>V</u> iew	<u>M</u> anage 🕜		A 🗖 🖥 🖉		
<mark> ∧ Pro<u>c</u>. Spectr</mark>	um 🚽 🐴 Adjust Phase 🚽	🔥 Calib. A <u>x</u> is 👻 🖇	<mark>¦%</mark> Pick P <u>e</u> aks ↓ ∫ <u>I</u> nte	grate 🗢 🛛 A <u>d</u> vanced 🗢			
Compute spectrum 1-click operation! options, click on th choose Configure	n from raw data (proc1d y). To change the processing e down arrow at right and Standard Processing.		·····································	TS EFP			
Browser Last50 Groups	Acquisition finished: /home/sb	ai/mynmr/demo/2/pdata/1			r 🛛 🖂		
የ-⊡ mynmr	Spectrum ProcPars AcquP	ars Title PulseProg	Peaks Integrals Sample	Structure Plot Fid Acqu			
<pre>> □ CCB</pre>	demo						
	200	150	100	50	0 [ppm]		
					s (ppm)		
Acquisition inform	tion Fid Flash Lock	Sample Shim Coil Temperature 303 K	POWCHK Sample Temperat	ure Spooler BSMS statu queued: 0 delayed: 0 Autoshim © L	is message Time 14:28:17 Nov 08		

Figure 32: Processed ¹³C Spectrum

6. Repeat steps 16 and 17 in section 3 (on page 14) to perform peak picking before printing the spectrum as described in step 21 of the section 3 (on page 16).

- In the BSMS panel, click on green button "Spin" to stop sample spinning, click on Lock ON-OFF button to turn off frequency lock, and finally click on green button "Lift" to eject sample. After the sample is removed from the magnet, click on Lift again to turn off the lifting air.
- 8. Exit the program as described in step 23 of section 3.