

The Liu Laboratory protocol — PSII-core Spinach  
Arts & Sci Washington University in St. Louis

2-2003

Oxygen-Evolving Reaction Center Core Preparation  
(Modified from Ghanotakis et. al., 1987)

Chl assay first.

1. Thaw out frozen OEP on ice, spin 35,000 x g for 30 min in the SS-34 head. Resuspend pellet to 2.5 mg/ml chl in SMN (or use fresh OEP finally suspended in SMN to 2.5 mg/ml chl). Measure volume. Also thaw out the stock 10%  $\beta$ -D-dodecyl maltoside (DM).

2. Measure an equal volume of Solution B - OGP. Add solid octyl glucopyranoside (OGP) to 2.06% (w/v), dissolve on rocker. Cool Solution B + OGP and add to OEP suspension. Incubate 10 min on ice.  $\checkmark$   $\sim 70\text{mM}$

in large flask (to accommodate increasing volume).  
3. To the total volume of solubilized membranes add two volumes of Solution C. Incubate 5 min, spin at 35,000 x g for 90 min in the SS-34 head. This will precipitate the LHC from the solution. Collect the supernatant, this contains the reaction center complex.  $20\text{ml} \times$

4. To the total volume of crude reaction center complex, add two volumes of Solution D and incubate for 5 min. Spin at 25,500 x g for 60 min in the GSA rotor. Carefully remove the bottles and CAREFULLY decant, and discard, the supernatant. The reaction center complex is found in the pellet. Remove as much supernatant as possible with a Pasture pipette. Resuspend the pellets in a small amount of Solution E (1-2 ml). Measure chl and bring to 1.0 mg/ml chl. Measure the volume using an ice cold 10 ml graduated cylinder.  $\text{stop here for}$

5. Transfer to a cold SS-34 centrifuge tube and add DM from a 10% stock solution to a final concentration of 0.5%. Incubate 60 min. Spin for 5 min at 35,000 x g in the SS-34 rotor. Carefully remove the supernatant and load onto a Sepharose 6B column previously equilibrated with 10 mM Mes-NaOH, pH 6.0, 5 mM NaCl, 5 mM CaCl<sub>2</sub>, 0.025% DM. Elute column at 1.5 ml/min.  $0.025$

6. As the green band approaches the bottom of the column, connect the fraction collector and collect 3.33 min fractions (5 ml). After the green band has eluted measure the absorbance of the fractions at 674 nm after diluting 100  $\mu$ l of fraction with 1 ml of elution buffer. Plot data as absorbance vs. volume. Discard the leading shoulder and pool the main band.

7. Fully concentrate in centricon 100. This is the Reaction Center Core Complex.

$\frac{2.06}{100} \cdot 3$   
 $\frac{2}{100} \cdot 3$   
 $\frac{6}{100}$   
 $\frac{6}{100}$   
 $\frac{6}{100}$   
6g  
6g

**The Liu Laboratory protocol — PSII-core Spinach  
Arts & Sci Washington University in St. Louis**

**Solutions**

1. SMN

2. Solution B - OGP

1 M sucrose

50 mM Mes-NaOH, pH 6.0

0.8 M NaCl

10 mM CaCl<sub>2</sub>

3. Solution C

1 M sucrose

50 mM Mes-NaOH, pH 6.0

0.4 M NaCl

5 mM CaCl<sub>2</sub>

4. Solution D

50 mM Mes-NaOH, pH 6.0

10 mM NaCl

5 mM CaCl<sub>2</sub>

5. Solution E

0.4 M sucrose

50 mM Mes-NaOH, pH 6.0

15 mM NaCl

5 mM CaCl<sub>2</sub>

6. Column Elution Buffer

10 mM Mes-NaOH, pH 6.0

5 mM NaCl

5 mM CaCl<sub>2</sub>

0.025% DM

SMN	400 mM Sucrose
	50 mM Mes-NaOH 6.0
	15 mM NaCl
	$\frac{0.15 \times 0.1}{0.15 + 0.1} = 0.5\text{M}$

$$\text{Solut. #2} \quad \text{Vol needed} \quad 0.1 \times 1 \times 342 = 34.2$$

$$\text{Solut. #3} \quad 200 \text{ uL} \quad 5 \times x = 0.8 \cdot 0.1 \\ = \frac{0.08}{0.08} = 1$$

$$\text{Solut. #4} \quad 600 \text{ uL} \quad 1 \cdot x = 0.1 \cdot 0.01 \\ x = 0.001$$

0.2 + 0.2

$$0.4 + 0.2 = 5 \cdot x \\ \frac{0.06}{0.06} = 1$$

$$0.6 \cdot 0.01 = 5 \cdot x \\ \frac{0.006}{0.006} = 1$$

$$0.6 \cdot 0.005 = 1 \cdot x \\ 0.0030 = 1$$

$$0.05 + 0.4 = 0.02$$

$$0.005 + 0.05 = 1 \cdot x$$

$$0.0055 = 1$$

Solution B:

100 ml
Sucrose 34.2 g
0.5 MES 6.0 10 ml
5 M NaCl 16 ml
1 M CaCl <sub>2</sub> 1 ml

Solution C:

200 ml
Sucrose 68.4 g
0.5 MES 6.0 20 ml
5 M NaCl 16 ml
1 M CaCl <sub>2</sub> 1 ml

Solution D:

600 ml
0.5 MES 60 ml
5 M NaCl 1.2 ml
1 M CaCl <sub>2</sub> 3 ml

Solution E:

50 ml
Sucrose 6.8 g
0.5 MES 5 ml
5 M NaCl 15 ml
1 M CaCl <sub>2</sub> 2.5 ml