

Oxidation of Alcohols

Chemicals and Equipment Needed

- 1-propanol – **R3**
- 2-propanol (isopropanol) – **R3**
- *t*-butanol (2-methyl-2-propanol) – **R1**
- Acidic $\text{K}_2\text{Cr}_2\text{O}_7$ solution – **K3**
 - Recipe follows
- 3 small crystalizing dishes – **P3**
- Plastic petri dishes – **P3**
- 3-10 mL graduated cylinders – **Q3**
- Stoppers for the cylinders (#00)– **U3**
- 3 glass rods – **U1**
- Plexiglass square – **J**
- Transparency from this packet

Hazards

- Alcohols are flammable and volatile solutions. Take care when handling, and pour them in the hood.
- The acidic $\text{K}_2\text{Cr}_2\text{O}_7$ solution is ~2 M H_2SO_4 . Wear gloves when handling.

Preparation

- Fill one cylinder with 3 mL 1-propanol, label **1°**. Fill the second with 3 mL 2-propanol, label **2°**. Fill the third with 3 mL *t*-butanol, label **3°**. Stopper all cylinders, and parafilm.
- Add 30 mL $\text{K}_2\text{Cr}_2\text{O}_7$ solution to each crystalizing dish, lid and label.
- On delivery, place the transparency and plexiglass square on the overhead, and the dishes on top (leave the H_2O space blank). Set the cylinders and glass rods next to the overhead.

Presentation

- Add the 1° and 2° alcohols to their respective dishes, then add the 3° alcohol to the proper dish.
 - You may wish to have a student time the reactions.
 - You should see a color change in the 1° and 2° dishes, but the total reaction time is 80 s for the 1° dish and 120 s for the 2° dish. The 3° alcohol does not react.

Discussion

- The 1° alcohol is oxidized to the aldehyde, then to the acid:
1-propanol → propanal → propionic acid
[ethanol → ethanol (acetaldehyde) → ethanoic acid (acetic acid)]
- The 2° alcohol is oxidized to the ketone:
2-propanol → propanone (acetone)
- As the alcohol is oxidized, the dichromate is reduced to Cr(III) ion, which will appear green, blue, or violet in these solutions:
$$\text{Cr}_2\text{O}_7^{2-}(\text{aq}) + 14 \text{H}^+(\text{aq}) + 6 \text{e}^- \rightarrow 2 \text{Cr}^{3+}(\text{aq}) + 7 \text{H}_2\text{O}(\text{g})$$
- The 3° alcohol cannot be oxidized, so it does not react:
t-butanol → no reaction
- NB: the dichromate oxidation of ethanol is the basis of the breathalyzer test.

Clean-Up

- Combine solutions in a beaker, and neutralize with 6 M NaOH. Pour into white waste container, along with any rinses.

NOTES:

Acidic $\text{K}_2\text{Cr}_2\text{O}_7$ solution: Mix 60 mL concentrated H_2SO_4 with 2.0 g $\text{K}_2\text{Cr}_2\text{O}_7$ (O2), dilute to 500 mL with d- H_2O . Final solution is ~2 M H_2SO_4

Lecturer Instructions

Presentation

- Add the 1° and 2° alcohols to their respective dishes, then add the 3° alcohol to the proper dish.
 - You may wish to have a student time the reactions.
 - You should see a color change in the 1° and 2° dishes, but the total reaction time is 80 s for the 1° dish and 120 s for the 2° dish. The 3° alcohol does not react.

Discussion

- The 1° alcohol is oxidized to the aldehyde, then to the acid:
1-propanol → propanal → propionic acid
[ethanol → ethanol (acetaldehyde) → ethanoic acid (acetic acid)]
- The 2° alcohol is oxidized to the ketone:
2-propanol → propanone (acetone)
- As the alcohol is oxidized, the dichromate is reduced to Cr(III) ion, which will appear green, blue, or violet in these solutions:
$$\text{Cr}_2\text{O}_7^{2-} (\text{aq}) + 14 \text{H}^+ (\text{aq}) + 6 \text{e}^- \rightarrow 2 \text{Cr}^{3+} (\text{aq}) + 7 \text{H}_2\text{O} (\text{g})$$
- The 3° alcohol cannot be oxidized, so it does not react:
t-butanol → no reaction
- NB: the dichromate oxidation of ethanol is the basis of the breathalyzer test.