

## LCM Preparation for Frozen Sections

### 1. Sample Freezing in OCT

- 1.1. Prepare your 2-methylbutane solution.
- 1.2. Pour 2-methylbutane in a metal bowl or plastic beaker. Add small pieces of dry ice into the solution to cool the solution down to  $-40^{\circ}\text{C}$  to  $-50^{\circ}\text{C}$ .
- 1.3. Collect dissected tissue specimen. Place tissue on paper towel to get rid of excess liquid before you put the tissue into the O.C.T.
- 1.4. Choose proper size of cryomold for tissue specimen.
- 1.5. Add a thin layer of O.C.T. compound (Optimal Cutting Temperature) to cryomold and place tissue specimen in desired orientation. Completely fill the mold with OCT. **There should be no bubbles and the tissue should be completely immersed in OCT.**
- 1.6. Hold the mold with hemostat and let the bottom of the mold touch the cold 2-methylbutane.
- 1.7. Place cryomold in cold ( $-40^{\circ}\text{C}$ ) 2-methylbutane solution when it is completely white and let it stay in the cold bath for another 5 minutes.
- 1.8. Wrap the frozen block with labeled foil. Store frozen sample in the cryomold in a  $-80^{\circ}\text{C}$  freezer or temporarily on dry ice.
- 1.9. Repeat these steps if you have more than one sample.

### 2. Slide Preparation

- 2.1. Pre-cool the cryostat to the temperature recommended by the manufacturer for the specimen you are preparing ( $-20^{\circ}\text{C}$  for most of tissues or  $-30^{\circ}\text{C}$  to  $-35^{\circ}\text{C}$  for lipid rich tissues).
- 2.2. Remove and discard old microtome blade. Wipe down the knife holder and anti-roll plate in the cryostat with 100% ethanol to avoid sample cross-contamination.
- 2.3. Install a new disposable microtome blade in the cryostat.
- 2.4. Set cutting thickness to desired thickness (**4-10  $\mu\text{m}$** ). Set cutting angle to **10** degrees (HM505EVP).
- 2.5. Transfer the cryomold containing the specimen from the  $-80^{\circ}\text{C}$  freezer to the cryostat, transporting on dry ice if necessary.
- 2.6. Wait a minimum of **10 minutes** for the specimen to equilibrate with the temperature of the cryostat.
- 2.7. Mount specimen to the specimen holder (chuck) with OCT. After frozen specimen trimming cut **4-10  $\mu\text{m}$**  sections.
- 2.8. Sections should be mounted within the membrane of the PEN-Membrane slides. Place slide immediately into slide box on dry ice (or  $-20^{\circ}\text{C}$  or  $-80^{\circ}\text{C}$  freezer). Do not allow slide to dry at room temperature.

- 2.9. If cutting more than one specimen, use a new disposable microtome blade for each one. In addition, wipe down knife holder and anti-roll plate with **100% ethanol** in between each specimen to avoid cross-contamination.
- 2.10. Proceed immediately to the “Staining and Dehydration” segment of the protocol or store at **-80°C** for up to two months.

***Note:** Frequent cycling of the tissue block from -80° C to -20° C for cryosectioning may accelerate RNA degradation. For best results, cut and mount a sufficient number of sections for two months' use during one cryosectioning session. Store the mounted sections at -80° C until needed.*

### 3. Alternative Protocol for H&E Staining and Dehydration

***Note:** Carry out the “H&E Staining and Dehydration” protocol with only one slide at a time. Change all solutions in the plastic slide jars between each case (maximum of 5 slides) to avoid cross contamination. Do not reuse solutions. Do not transfer solutions back into their original bottles. If you plan to reuse jars, discard all water-based solutions upon completion of staining. Clean them with 100% ethanol, followed by distilled water, RNase AWAY or RNase ZAP, then distilled water and allow drying completely in the hood).*

- 3.1. Label and fill seven slide coplin jars as follows:

	Coplin jar (glass or plastic)	Time
<b>a</b>	<b>70% Ethanol (made with DEPC water)</b>	<b>minimum of 30 seconds @ room temp (-20°C for long term storage)</b>
<b>b</b>	<b>DEPC water</b>	<b>30 seconds</b>
<b>c</b>	<b>Mayer's Hematoxylin</b>	<b>2 minutes</b>
	<b>Gill II or III</b>	<b>10-30 seconds</b>
	<b>Methylene Blue (better for DNA)</b>	<b>2 minutes</b>
<b>d</b>	<b>DEPC water</b>	<b>30 seconds</b>
<b>e</b>	<b>Alcoholic Eosin</b>	<b>10 seconds</b>
<b>f</b>	<b>95% Ethanol</b>	<b>30 seconds</b>
<b>g</b>	<b>100% Ethanol</b>	<b>30 seconds</b>
<b>h</b>	<b>Air dry</b>	<b>1-5 minutes</b>

- 3.2. Fill the labeled coplin slide jars with **50 ml** of the appropriate solution (see table).
- 3.3. Remove one slide at a time from the slide box on dry ice (or -20°C or -80°C freezer), and place it in slide jar “a” containing **70% ethanol** for **30 seconds**. Use cover glass forceps to transfer slides from jar to jar.
- 3.4. Transfer the slide to jar “b” containing **DEPC water** for **30 seconds**.
- 3.5. Transfer the slide to jar “c” containing hematoxylin solution for appropriate time (**see table**).
- 3.6. Place the slide in jar “d” containing **DEPC water** for **30 seconds**.
- 3.7. Transfer the slide to jar “e” containing **alcoholic eosin** solution for **10 seconds**.
- 3.8. Transfer the slide to jar “f” containing **95% ethanol** for **30 seconds**.
- 3.9. Transfer the slide to jar “g” containing **100% ethanol** for **30 seconds**.

- 3.10. Air-dry slide for Leica AS LMD for **1-2 minutes** and start laser microdissection immediately.
- 3.11. Discard all used staining and dehydration solutions according to your procedures.