

# Microvessel Analysis Algorithm User's Guide



# Microvessel Analysis Algorithm User's Guide

This document applies to eSlide Manager Release 12.3 and later.

## Copyright Notice

- ▶ Copyright © 2008-2015 Aperio. All rights reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. Aperio is a registered trademark of Leica Biosystems Imaging, Inc. in the USA and other countries.

## Customer Resources

- ▶ For the latest information on Leica Biosystems Aperio ePathology products and services, please visit [www.LeicaBiosystems.com/ePathology](http://www.LeicaBiosystems.com/ePathology).

## Disclaimers

- ▶ Use normal care in maintaining and using Aperio ePathology servers. Interrupting network connections or turning off the servers while they are processing data (such as when they are analyzing eSlides or generating an audit report) can result in data loss.
- ▶ This manual is not a substitute for the detailed operator training provided by Leica Biosystems Imaging or for other advanced instruction. Leica Biosystems Imaging Field Representatives should be contacted immediately for assistance in the event of any instrument malfunction. Installation of hardware should only be performed by a certified Leica Biosystems Imaging Service Engineer.
- ▶ ImageServer is intended for use with eSlides created by scanning glass slides with the scanner. Educators will use Aperio ePathology software to view and modify eSlides in Composite WebSlide (CWS) format.

## Patents

- ▶ Aperio ePathology products are protected by U.S. Patents: 6,711,283; 6,917,696; 7,035,478; 7,116,440; 7,257,268; 7,428,324; 7,457,446; 7,463,761; 7,502,519; 7,518,652; 7,602,524; 7,646,496; 7,738,688 and licensed under one or more of the following U.S. Patents: 6,101,265; 6,272,235; 6,522,774; 6,775,402; 6,396,941; 6,674,881; 6,226,392; 6,404,906; 6,674,884; and 6,466,690.

## Contact Information – Leica Biosystems Imaging, Inc.

Headquarters	Customer Support	General Information
 Leica Biosystems Imaging, Inc. 1360 Park Center Drive Vista, CA 92081 USA  Tel: +1 (866) 478-4111 (toll free) Direct International Tel: +1 (760) 539-1100	US/Canada Tel: +1 (866) 478-3999 (toll free) Direct International Tel: +1 (760) 539-1150  US/Canada/Worldwide Email: TechServices@LeicaBiosystems.com	US/Canada Tel: +1 (866) 478-4111 (toll free) Direct International Tel: +1 (760) 539-1100  Email: ePathology@LeicaBiosystems.com

# Customer Service Contacts

Please contact the office for your country for technical assistance.

**Australia:**

96 Ricketts Road  
Mount Waverly, VIC 3149  
AUSTRALIA  
Tel: 1800 625 286 (toll free)  
Between 8:30 AM-5 PM, Monday-Friday, AEST  
Email: lbs-anz-service@leicabiosystems.com

**Austria:**

Leica Biosystems Nussloch GmbH  
Technical Assistance Center  
Heidelberger Strasse 17  
Nussloch 69226  
GERMANY  
Tel: 0080052700527 (toll free)  
In-country Tel: +43 1 486 80 50 50  
Email: support.at@leicabiosystems.com

**België/Belgique:**

Tel: 0080052700527 (toll free)  
In-country Tel: +32 2 790 98 50  
Email: support.be@leicabiosystems.com

**Canada:**

Tel: +1 866 478- 999 (toll free)  
Direct International Tel: +1 760 539 1150  
Email: TechServices@leicabiosystems.com

**China:**

17F, SML Center No. 610 Xu Jia Hui Road, Huangpu District  
Shanghai, PRC PC:200025  
CHINA  
Tel: +86 4008208932  
Fax: +86 21 6384 1389  
Email: service.cn@leica-microsystems.com  
Remote Care email: tac.cn@leica-microsystems.com

**Danmark:**

Tel: 0080052700527 (toll free)  
In-country Tel: +45 44 54 01 01  
Email: support.dk@leicabiosystems.com

**Deutschland:**

Leica Biosystems Nussloch GmbH  
Technical Assistance Center  
Heidelberger Strasse 17  
Nussloch 69226  
GERMANY  
Tel: 0080052700527 (toll free)  
In-country Tel: +49 6441 29 4555  
Email: support.de@leicabiosystems.com

**Eire:**

Tel: 0080052700527 (toll free)  
In-country Tel: +44 1908 577 650  
Email: support.ie@leicabiosystems.com

**España:**

Tel: 0080052700527 (toll free)  
In-country Tel: +34 902 119 094  
Email: support.spain@leicabiosystems.com

**France:**

Tel: 0080052700527 (toll free)  
In-country Tel: +33 811 000 664  
Email: support.fr@leicabiosystems.com

**Italia:**

Tel: 0080052700527 (toll free)  
In-country Tel: +39 0257 486 509  
Email: support.italy@leicabiosystems.com

**Japan:**

1-29-9 Takadannobaba, Sinjuku-ku  
Tokyo 169-0075  
JAPAN

**Nederland:**

Tel: 0080052700527 (toll free)  
In-country Tel: +31 70 413 21 00  
Email: support.nl@leicabiosystems.com

**New Zealand:**

96 Ricketts Road  
Mount Waverly, VIC 3149  
AUSTRALIA  
Tel: 0800 400 589 (toll free)  
Between 8:30 AM-5 PM, Monday-Friday, AEST  
Email: lbs-anz-service@leicabiosystems.com

**Portugal:**

Tel: 0080052700527 (toll free)  
In-country Tel: +35 1 21 388 9112  
Email: support.pt@leicabiosystems.com

**Sweden:**

Tel: 0080052700527 (toll free)  
In-country Tel: +46 8 625 45 45  
Email: support.se@leicabiosystems.com

**Switzerland:**

Tel: 0080052700527 (toll free)  
In-country Tel: +41 71 726 3434  
Email: support.ch@leicabiosystems.com

**United Kingdom:**

Tel: 0080052700527 (toll free)  
In-country Tel: +44 1908 577 650  
Email: support.uk@leicabiosystems.com

**USA:**

Tel: +1 866 478 3999 (toll free)  
Direct International Tel: +1 760 539 1150  
Email: TechServices@leicabiosystems.com

# Contents

- 1 Introduction ..... 7**
  - About This Guide.....7
  - About the Microvessel Analysis Algorithm .....7
  - Prerequisites .....8
  - Intended Use.....8
  - Installing the Algorithm .....8
  - For More Information .....8
  - Acknowledgments .....9
  
- 2 Tuning Parameters..... 10**
  - Tuning Algorithm Parameters .....10
  - Algorithm Parameters.....10
  - Algorithm Results.....12
  
- 3 Microvessel Analysis Process ..... 14**
  - Microvessel Analysis Definitions.....14
  - Tuning the Microvessel Analysis Algorithm Macro .....15
  
- 4 Color Deconvolution..... 16**
  - Setting the Stain Color Vectors .....17
  - Checking the Deconvolution.....18
  
- 5 Defining Light and Dark Staining Thresholds ..... 19**
  - Vessel Thresholding .....19
  
- 6 Completing Vessels ..... 21**
  - Region Joining Parameter .....21
  - Vessel Completion Parameter.....22
  - Completing the Vessel .....22
  
- 7 Analysis Metrics..... 25**
  - Selecting Maximum/Minimum Vessel Size to Analyze.....25
  - Additional Analysis Modes .....26
  - Generating Histograms .....26

Saving and Using the Algorithm Macro.....	28
Example - Oncology Research with Xenografts .....	28
For More Information .....	28
<b>Index .....</b>	<b>29</b>
<b>Symbols.....</b>	<b>31</b>

# 1

## Introduction

This chapter introduces the microvessel analysis algorithm. For general information on using an algorithm, please see the *Aperio Image Analysis User's Guide*.



*The primary source for information on creating, testing, and saving image analysis macros is the **Aperio Image Analysis User's Guide**. That guide also contains details on running image analyses from within eSlide viewers and batch analyses from within eSlide Manager, and viewing and exporting analysis results.*



*The Aperio Image Analysis Workstation provides a streamlined image analysis workflow on your local workstation; if you are using that product, please refer to the **Aperio Image Analysis Workstation User's Guide** for instructions on creating, testing, and saving image analysis macros. That guide also contains details on running image analyses from within ImageScope, running batch analyses, and viewing and exporting analysis results.*

### About This Guide

This guide for this image analysis algorithm discusses how to set the algorithm parameters to suit your image analysis needs. After tuning the parameters, you will save the settings as an algorithm macro. The macro can then be used by you and other users to analyze specific eSlides (digital slide images).

The guide works in concert with the *Aperio Image Analysis User's Guide* or the *Aperio Image Analysis Workstation User's Guide* to present the complete picture of Aperio image analysis.

### About the Microvessel Analysis Algorithm

Analyzing capillary and microvessel infrastructure is critical in many disease areas. From in vivo angiogenesis assays in cardiovascular research to vascular morphology changes in oncology xenografts to microcapillary measurements in rodent models in ophthalmology, researchers need to understand and quantify vascular architecture. With whole slide scanning and the microvessel analysis algorithm, a tedious and subjective manual counting task can be automated and standardized.

The microvessel analysis algorithm detects and quantifies microvessels on slides stained with endothelial markers (for example, CD31, CD34, and Factor VIII).

Microvessel analysis typically starts with a Factor VIII or CD31 stained slide, scanned at 20x or 40x. A dual-level thresholding process screens out non-specific staining. Regions are then assembled into vessels, and individual statistics are calculated on each vessel. You can choose which statistic to graph, and then can generate a summary report.

## Prerequisites

Because Aperio eSlides are by design high resolution and information rich, for best results you should use a high quality monitor to view them. Make sure the monitor is at the proper viewing height and in a room with appropriate lighting. We recommend any high quality LCD monitor meeting the requirements recommended in the *Aperio ePathology System Requirements*.

## Intended Use

For research use only. Not for use in diagnostic procedures.

Algorithms are intended to be used by trained pathologists who have an understanding of the conditions they are testing for in running the algorithm analysis.

Each algorithm has input parameters that must be adjusted by an expert user who understands the goal of running the analysis and can evaluate the algorithm performance in meeting that goal.

You will adjust (tune) the parameters until the algorithm results are sufficiently accurate for the purpose for which you intend to use the algorithm. You will want to test the algorithm on a variety of images so its performance can be evaluated across the full spectrum of expected imaging conditions. To be successful, it is usually necessary to limit the field of application to a particular tissue type and a specific histological preparation. A more narrowly defined application and consistency in slide preparation generally equates to a higher probability of success in obtaining satisfactory algorithm results.

If you get algorithm analysis results that are not what you expected, please see the “Troubleshooting” appendix in the *Aperio Image Analysis User’s Guide* for assistance.

## Installing the Algorithm

In most cases you install the algorithm on the eSlide Manager server only. (In fact, Technical Services may install the algorithm for you on your server.) This is because you typically fine-tune and save the algorithm parameters on the eSlide Manager server. For the rare case that you need to fine-tune the algorithm parameters on your local workstation or when using a local image, refer to the *Aperio Image Analysis User’s Guide* for installation and use instructions.

## For More Information

For a quick reference to the microvessel analysis algorithm parameters, refer to “Chapter 2: Tuning Parameters” on page 10.

For details on using the microvessel analysis algorithm, begin with “Chapter 3: Microvessel Analysis Process” on page 14. The remaining chapters give details on using the algorithm.

See the *Aperio Image Analysis User’s Guide* for information on:

- ▶ Installing an algorithm
- ▶ Creating a new algorithm macro or modifying an existing one
- ▶ Saving or exporting a macro
- ▶ Selecting the areas of an eSlide to analyze



- ▶ Running an analysis on a single eSlide through the eSlide viewer or running an analysis on one or more eSlides using eSlide Manager batch analysis
- ▶ Viewing analysis results quantitatively and visually
- ▶ Exporting analysis results

For details on using ImageScope to view eSlides, see the *ImageScope User's Guide*.

If the analysis results are not what you expect, see the *Aperio Image Analysis User's Guide* section on troubleshooting for assistance.

## Acknowledgments

Leica Biosystems Imaging thanks the following researchers for their suggestions and for donating images used in this guide:

- ▶ Mouse skin containing a subcutaneous matrigel plug immunostained with Factor VIII-related antigen provided by Dr. Reginald Valdez, DVM, MS, PhD, DACVP, Novartis Institutes for BioMedical Research, Inc.
- ▶ James Deeds, Novartis NIBR for many helpful suggestions.
- ▶ Dr. Maher Albitar, MD, Ph.D. Quest Diagnostics Nichols Institute.
- ▶ Dr. Judah Folkman, MD for taking the time to discuss the algorithm, his suggestions on the vessel perimeter measurements for a correlative measure of oxygenation, and for his clear review articles in the field that concisely describe what is required for improved future monitoring of angiogenesis.

# 2

## Tuning Parameters

This chapter contains information on all microvessel analysis algorithm inputs and outputs. See the following chapters for details on using the algorithm.

### Tuning Algorithm Parameters

To create a macro for the algorithm or to modify an existing macro:


1. In ImageScope, open an eSlide you want to use to tune the algorithm parameters. (Refer to the *ImageScope User's Guide* for instructions.)
2. In ImageScope, open the algorithm and choose to create a new macro or open an existing macro to modify it. (Refer to the *Aperio Image Analysis User's Guide* for instructions.)

You now see the parameters for the algorithm (these are listed later in this chapter).

3. Now adjust the parameters in the Analysis window as discussed below and move the Algorithm Tuning window on the image to see a mark-up image of the results.

After adjusting the input parameters and selecting the output parameters as discussed in the following sections, move the Algorithm Tuning window to various areas on the eSlide to see an approximation of the analysis results. The results appear in the Annotations window and as a mark-up image in the ImageScope main window.

You can open the ImageScope Annotations window to see the algorithm tuning results displayed numerically.

When saving the algorithm macro by clicking , you can choose whether to save the macro locally on your workstation or, if connected to eSlide Manager, to save the macro remotely on eSlide Manager. You are asked to supply a name for the macro that will help you identify it in the future. If saving the macro to eSlide Manager, you are also asked to specify which data group the macro will be associated with; only eSlide Manager users who have permission to use that data group will be able to choose that macro for analysis.

### Algorithm Parameters

Microvessel analysis algorithm performance is controlled by a set of input parameters, which determine many different types of analysis.

- ▶ **View Width** – Width of processing box.
- ▶ **View Height** – Height of processing box.
- ▶ **Overlap Size** – Size of the overlap region for each view. This should be at least as big as the largest vessel of interest.
- ▶ **Image Zoom** – Zoom level to be used; a higher zoom results in faster algorithm run but less accurate results.

- ▶ **Markup Compression Type** – This sets the compression type for the algorithm mark-up image. Choose better compression if you need the image for a special purpose.
- ▶ **Compression Quality** – A higher quality takes longer and yields larger files. This selection does not apply to all compression types.
- ▶ **Markup Image Type** – There are two types of markup Images – Tuning and Analysis Markups. Tuning markups help in fine-tuning various input parameters before the vessel analysis. Only a few numeric results are valid for these markups and the rest of the results are zero. Analysis markups show different results of the microvessel analysis algorithm. All the numeric results are valid and identical for the analysis markups.
- ▶ **Mode** – Choose whether to include/exclude incomplete vessels and other stained objects in the numeric results.
- ▶ **Filtering / Smoothing Level** – Specifies amount of filtering on the deconvolved vessel channel. Ranges from 0 – 10 with zero resulting in no filtering and 10 the maximum level of filtering.
- ▶ **Dark Staining Threshold** – First of the two parameters in segmenting the deconvolved vessel channel. (0 -255). Vessels don't stain uniformly at all locations depending on the "cut." This intensity corresponds to the darker or more intense regions of the vessels
- ▶ **Light Staining Threshold** – Second parameter in segmenting the deconvolved vessel channel. This threshold corresponds to the majority of vessel regions that stain at least as light or intense as this value.
- ▶ **Region Joining Parameter** – Distance in microns between pieces of vessel walls that may be joined into a complete vessel.
- ▶ **Vessel Completion Parameter** – Distance in microns between the ends of a vessel wall that may be completed into a vessel.
- ▶ **Minimum Vessel Area Threshold** – Any vessel which has an internal area (in  $\mu\text{m}^2$ ) less than this value is eliminated from numeric results.
- ▶ **Maximum Vessel Area Threshold** – Any vessel which has an internal area (in  $\mu\text{m}^2$ ) more than this value is eliminated from numeric results.
- ▶ **Maximum Vessel Wall Thickness** – Vessel Walls might seem much thicker than their actual size depending on the way they are cut. A non-orthogonal cut usually result in thicker looking Vessel Walls. This parameter restricts the maximum possible vessel wall thickness.
- ▶ **Output Histogram** – Choice of output statistic for histogram analysis. "None" option can be chosen for no histogram analysis.
- ▶ **Histogram Start Value** – Specifies the starting point for histogram analysis. Any output statistic value less than the start value (in either  $\mu\text{m}$  or  $\mu\text{m}^2$ ) is not included in the histogram results.
- ▶ **Histogram End Value** – Specifies the end point for histogram analysis. Any output statistic value greater than the end value (in either  $\mu\text{m}$  or  $\mu\text{m}^2$ ) is not included in the histogram results. Value of "-1" forces the histogram analysis on all output statistic values.
- ▶ **Number of Bins** – Specifies the number of bins in histogram analysis. All the histograms are uniformly spaced.
- ▶ **Endothelial Stain Red, Green, Blue Components** – Normalized optical density values for Endothelial Stain. Default for endothelial stain is DAB. These values are calculated by a calibration procedure described earlier in this guide.

- ▶ **Background Stain (1) Red, Green, Blue Components** – Normalized optical density values for Background channel (1). Default for background channel (1) is Hematoxylin. These values are calculated by a calibration procedure described earlier in this guide.
- ▶ **Background Stain (2) Red, Green, Blue Components** – Normalized optical density values for Background channel (2). Default for background channel (2) is “None / No Stain”. These values are calculated by a calibration procedure described earlier in this guide.
- ▶ **Clear Area Intensity** – This is the intensity for a clear area on the slide. This value is always 240 for images generated by an Aperio scanner.

## Algorithm Results

Many of the algorithm results are discussed in the following chapters.

- ▶ **Number of Vessels** – Total number of vessels detected after microvessel analysis.
- ▶ **Total Analysis Area ( $\mu\text{m}^2$ )** – Total area that is analyzed in square-microns.
- ▶ **Total Stain Area ( $\mu\text{m}^2$ )** – Total area that has significant endothelial staining. These regions are segmented by the Dark and Light thresholds.
- ▶ **Average Stain Intensity** – Average stain intensity of all the endothelial stained pixels.
- ▶ **Microvessel Density ( $\mu\text{m}^2$ )** – Number of vessels per unit of analysis area in square-microns.
- ▶ **Mean Vessel Area ( $\mu\text{m}^2$ )** – Average of all the Vessel Areas in square-microns.
- ▶ **Median Vessel Area ( $\mu\text{m}^2$ )** – Median values of all the Vessel Areas in square-microns.
- ▶ **Standard Deviation of Vessel Area ( $\mu\text{m}^2$ )** – Standard deviation of all the Vessel Area in square-microns.
- ▶ **Mean, Median and Standard Deviations of Vessel Perimeter ( $\mu\text{m}$ ), Lumen Area ( $\mu\text{m}^2$ ), Vascular Area ( $\mu\text{m}^2$ ) and Vessel Wall Thickness ( $\mu\text{m}$ )** – Similar statistics for the parameters as above.
- ▶ **Average Red OD** – Average OD (optical density) of the Red component for all pixels analyzed. This value is used for calibration of input color specifications.
- ▶ **Average Green OD** – Average OD (optical density) of the Green component for all pixels analyzed. This value is used for calibration of input color specifications.
- ▶ **Average Blue OD** – Average OD (optical density) of the Blue component for all pixels analyzed. This value is used for calibration of input color specifications.

## Histogram Results

For an example of a histogram of results, see *“Generating Histograms” on page 26*.

- ▶ **Total Number of Vessels** – Total number of vessels detected by the algorithm
- ▶ **Number of vessels in Histogram Analysis** – Number of vessels that are included in the histogram analysis based on the histogram start and end values specified by input parameters.
- ▶ **Overall Minimum Value ( $\mu\text{m}^2$ )** – Overall minimum value of the statistic which is chosen for histogram analysis. This value is not affected by the histogram start and end values.
- ▶ **Overall Maximum Value ( $\mu\text{m}^2$ )** – Overall maximum value of the statistic which is chosen for histogram analysis.

This value is not affected by the histogram start and end values.

- ▶ **Histogram Bin Centers - Frequency** – Histogram Analysis results showing a list of bin centers and their corresponding frequencies.

# 3

## Microvessel Analysis Process

This chapter discusses the steps in microvessel analysis and gets you ready to start using the algorithm.

There are five steps to microvessel analysis:

**1. Scan the slide (create an eSlide)**

Use the Aperio scanner to scan a glass slide to create an eSlide. The eSlide file is automatically saved in eSlide Manager. See the user's guide for your Aperio scanner for details on scanning.

**2. Color deconvolution**

If you know that the stain used for the slide is DAB, the default color vectors are correct so you can skip this step; otherwise, you need to set the stain color vectors. See *"Chapter 4: Color Deconvolution"* on page 16.

**3. Light and dark staining thresholds**

You need to set the light and dark staining threshold values to maximize selecting vessels without selecting background tissue. See *"Defining Light and Dark Staining Thresholds"* on page 19.

**4. Completing vessels**

Next you need to join the incomplete vessels created in the previous step to create complete vessels. See *"Completing Vessels"* on page 21.

**5. Analysis metrics**

Now that you have identified complete vessels, you can use the algorithm to output histograms for such information as vessel area, vessel perimeter, vascular area, and so on. See *"Chapter 7: Analysis Metrics"* on page 25.

### Microvessel Analysis Definitions


The following terms are used in this user's guide. At various steps in using the microvessel algorithm, you will be able to see graphical representations of these elements.

Term	Definition
Endothelial cell	Epithelial cell that lines the circulatory vessel.
Region	The intermediate object the algorithm constructs by putting together endothelial cells.
Vessel	Connection of regions into one object, typically with a lumen
Lumen	Inside area of a vessel.
Vessel wall	Endothelial cells that form the wall of the vessel.

Term	Definition
Vascular area	Sum of all endothelial cells.
Vessel area	Vascular area/2 + lumen area

## Tuning the Microvessel Analysis Algorithm Macro

Once you have selected an eSlide and opened it in ImageScope (as discussed in the *Aperio Image Analysis User's Guide*), you see the eSlide displayed in the main ImageScope window. (The following is a digital image of mouse skin that contains a subcutaneous matrigel plug which has been immunostained with Factor VIII-related antigen.)



**Image of mouse skin containing a subcutaneous matrigel plug immunostained with Factor VIII-related antigen provided by Dr. Reginald Valdez, DVM, MS, PhD, DACVP, Novartis Institutes for BioMedical Research, Inc.**

For details on moving around the image, changing its magnification, and using the drawing tools to select areas of the image, see the *ImageScope User's Guide*.

# 4

## Color Deconvolution

Color Deconvolution defines the stain color vectors so that endothelial cells will be correctly identified by the algorithm. If the slide was stained with DAB, the default color vectors are correct and you can skip this step by going directly to the next chapter.

By defining the stain color vectors, you are identifying to the microvessel analysis algorithm which color identifies the endothelial cells that make up the vessel walls.

Deconvolution separates the contribution from the endothelial stain (for example, CD31 or Factor VIII) and uses only this color for analysis.


The default color vector settings are for DAB as the endothelial stain and hemotoxylin as the background stain 1 (background stain 2 set to zero; that is, assumed to be nonexistent).

Endothelial Stain - Red Compon...	0.268
Endothelial Stain - Green Comp...	0.57
Endothelial Stain - Blue Compon...	0.776
Background Stain 1 - Red Com...	0.68
Background Stain 1 - Green Co...	0.57
Background Stain 1 - Blue Com...	0.44
Background Stain 2 - Red Com...	0
Background Stain 2 - Green Co...	0
Background Stain 2 - Blue Com...	0

These numbers must be changed if different stains are used. The color for each stain is calibrated separately, using a separate image for each stain having only that color present.

Look for several areas of the eSlide that are mostly stained with endothelial stain and select them by using the ImageScope drawing tools. In other words, pick an area of light staining of only this color. Avoid selecting darker, over-stained areas.





For this example, we selected several light brown areas using the ImageScope drawing tools on the ImageScope toolbar. The tools you can use to select regions are:



Freehand pen – Use to draw a free-form area of interest.



Negative freehand pen – Use to draw an area to exclude from the analysis. Note that you can use this in combination with the other drawing tools to first select an area of interest and then exclude areas within the selected area that you do not want to analyze.




Rectangle tool – Draws a rectangular area. If you want to select a square, hold down the Shift key while drawing.

## Setting the Stain Color Vectors

Create a new Microvessel Analysis algorithm macro or test an existing one. (See the *Aperio Image Analysis User's Guide* for details on creating a new macro or modifying an existing one.)

On the Analysis window:

1. From the **Mark-up Image Type** drop-down list, select **Tuning – Vessel or Endothelial Stain**.
2. The Algorithm Tuning window shows a visual result of the analysis.



- To see the numeric results of the analysis, go to the **View** menu and select **Annotations**. The Annotations window displays the color vector information in the OD (optical density) area:

Average Red OD	0.426768
Average Green OD	0.58289
Average Blue OD	0.691454

- Type the values from the **Endothelial Stain – Red OD**, **Endothelial Stain – Green OD**, and **Endothelial Stain – Blue OD** lines into the corresponding areas of the Algorithm window:

Endothelial Stain - Red Component	0.268
Endothelial Stain - Green Component	0.57
Endothelial Stain - Blue Component	0.776

Now you will test this macro on other eSlides that contain the stains you need to set the color vectors for.

- Open an eSlide that is predominately stained with the stain you need to determine the color vector for.
- Perform the steps from the section above for the additional stain.

## Checking the Deconvolution

After you have set the color vectors, you can double-check your results by opening an eSlide that contains the tissue you are interested in, drawing a region around a vessel, run the algorithm again in tuning mode, and the vessel walls should be clearly seen.

Leave the Analysis window open for the steps in the next chapter.

# 5

## Defining Light and Dark Staining Thresholds

In order to distinguish between endothelial cells and background tissue, you need to define the light and dark staining thresholds.

Defining the light and dark staining thresholds will identify endothelial cells.

- ▶ **Dark Staining Threshold** – All vessels must contain some region that stains darker than this threshold.
- ▶ **Light Staining Threshold** – All vessels must contain large regions that stain darker than this threshold.

By setting the threshold values correctly, you avoid single-point staining (for example, false stain uptake by non-endothelial cells).

The darkest value equals zero, and the lightest value equals 255.

### Vessel Thresholding

With the eSlide open in ImageScope and the microvessel analysis algorithm parameters displayed on the Analysis window:

1. Select **Tuning – Vessel Thresholding** in the **Mark-up Image Type** drop-down list.
2. Set both threshold values to the same value (try 210).
3. Move the tuning window to the area of the eSlide you want to test. Here is a mark-up image when light and dark thresholds are set to 210. As you can see, quite a few non-vessel areas are selected.



4. Now set the dark threshold to a new value to something less than 190. Keep running the algorithm in this tuning mode and inspecting the markup image.

5. Keep reducing the dark threshold value to eliminate non-vessel material, but without eliminating vessels. Here is a mark-up image when light threshold is set to 210 and dark threshold is set to 170:



6. Adjust the light staining threshold value to ensure that nearly all endothelial cells are detected but a minimum of false stained cells are detected.

Your goal for this process is to strike a balance between selecting vessels but not selecting much background tissue. It is okay at this point to have a few false stained cells as the next step can eliminate them.

When you have the algorithm parameters set so that you are satisfied with the results, leave the macro open for the next chapter.

# 6

## Completing Vessels

In this step, you join incomplete vessels into complete vessels.

This step looks for lumen and connects endothelial cells around the lumen to identify a vessel.

Two things can be set:

- ▶ **Region Joining Parameter** – Connects broken pieces of vessel walls; all pieces within the distance specified will be joined to the same vessel.
- ▶ **Vessel Completion Parameter** – Connects incomplete and open vessel walls; incomplete vessel walls with less than this distance of an opening may be completed into one vessel.

You will modify the Microvessel Analysis algorithm macro you created in the earlier chapters.

### Region Joining Parameter


Here is an example of when the Region Joining Parameter can be used to connect cells into a vessel. Notice that this incomplete vessel consists of three regions separated by gaps:



The Region Joining Parameter specifies the maximum distance of the gaps that can be detected and filled in to join two regions when more than one gap exists. Once the regions have been joined so that a single object is created with just one gap, the Vessel Completion Parameter is used to complete the vessel.

## Vessel Completion Parameter

Here is an example of when the Vessel Completion Parameter can be used to close an incomplete vessel:




In this example, a complete object exists with a single gap.

## Completing the Vessel

First we will measure the gaps in the vessel:

1. Select the ruler tool on the ImageScope toolbar: .
2. Choose the worst-case incomplete vessel (the gaps are the largest) and drag the ruler across the gaps in the vessel to measure them. The measurements in microns display on the image:



3. On the Analysis window, go to the **Mark-up Image Type** drop-down list and select **Analysis – Vessels**:
4. Type a value for the largest gap measurement into **Vessel Completion Parameter** and the next largest gap measurement into **Region Joining Parameter**.
5. Move the Algorithm Tuning window to the area you have annotated. If the vessel has been completed, you will see an image similar to this, which shows the lumen filling the vessel:



If you supply Vessel Completion Parameter and Region Joining Parameter values that are too large, the algorithm will join multiple vessels together and will also take much longer to run.

6. To see numeric information from the analysis, go to the ImageScope **View** menu and select **Annotations**. The Annotations window displays with information on the different tuning and analysis operations we have performed:

The screenshot shows the 'Annotations' window with the following data:


Layer Attributes	
Number of Vessels	2
Total Analysis Area (um2)	9548.6
Total Stain Area (um2)	1421.6
Average Stain Intensity	180.369
Microvessel Density -- Number of vessels per unit area (um2)	2.09454e-004
Mean Vessel Area (um2)	1423.5
Median Vessel Area (um2)	2739.
Standard Deviation of Vessel Area (um2)	1860.4
Mean Vessel Perimeter (um)	191.5
Median Vessel Perimeter	336.

Layer Regions			
Region	Length (um)	Area (um2)	Text
1	392.4	9595	



- Click on the different annotation layers in the Layers pane to show numeric data for that step (and to show the mark-up image for that step in the ImageScope window). For example, to see the ruler measurements again, click on the annotation layer for those measurements (in the example above, we have renamed that layer to "Ruler Measurements"), and you see both the numeric results (the measurements in the Layer Regions pane) and the mark-up image for that step.



- Leave the Analysis window open for the steps in the next chapter.



# 7

## Analysis Metrics

Now that complete vessels have been identified, you can use the microvessel algorithm to analyze microvessel distribution to create histograms.

By using the microvessel analysis algorithm parameters on the Algorithms window, you can select many different types of analysis.


Please see *“Chapter 2: Tuning Parameters” on page 10* for details on each parameter.

By selecting the appropriate algorithm parameters and mark-up image types, you can fine-tune the microvessel selection process (for example, including or excluding microvessels to analyze based on minimum or maximum size), output histograms for vessel area, vessel perimeter, vessel wall thickness, vascular area, and lumen area.

### Selecting Maximum/Minimum Vessel Size to Analyze

To fine-tune microvessel detection by selecting the minimum and maximum size vessels to analyze:

1. Draw freehand shapes around maximum and minimum size vessels . You see the area of those selections in the Annotations window:



2. On the Analysis window, set the maximum and minimum vessel size parameters to reflect your selections. For example, to exclude the smaller vessel selected, enter a minimum size of greater than 370; to include vessels as large or larger than the larger vessel selected, enter a maximum size greater than 4492.
3. Now the algorithm runs in tuning mode, and the smaller vessel is excluded from analysis, but the larger one is included as shown by the mark-up images below:



### Additional Analysis Modes

The different **Mark-up Image Type** selections will give you many different ways to look at microvessels. For example, choosing **Analysis – Lumen and Vascular Cells** results in a mark-up image like this:



You can also use the **Mode** drop-down list to include/exclude incomplete vessels from analysis:

- 0 - Exclude Incomplete Vessels & Other Stained Regions
- 1 - Include Incomplete Vessels & Other Stained Regions

### Generating Histograms

The **Output Histogram** drop-down list on the Algorithm window gives several different choices for the data to output as a histogram:

- ▶ 0 – None
- ▶ 1 – Vessel Area
- ▶ 2 – Vessel Perimeter
- ▶ 3 – Lumen Area
- ▶ 4 – Vascular Area


▶ 5 – Vessel Wall Thickness

There are additional parameters that affect how the histogram data is output:

- ▶ Histogram Start Value ( $\mu\text{m}$  or  $\mu\text{m}^2$ )
- ▶ Histogram End Value ( $\mu\text{m}$  or  $\mu\text{m}^2$ )
- ▶ Number of bins


Note that you can adjust the start and end values for the histogram to capture just the data you want. A histogram end value of -1 captures all data. You might, however, want to set this to the value of the largest vessel for which you want data (for example, to fine-tune the range of vessels reported on).

When you have selected a histogram type and run the algorithm, the histogram data is shown in the Layer Attributes pane of the Annotation window.

Now you can use the spreadsheet export button on the Layer Attributes pane, , to export the data into an Excel spreadsheet where you can see the data in graphical form.

For example:

<b>Statistic</b>	<b>Tumor</b>	<b>Normal</b>
Number of Vessels	704	352
Total Analysis Area ( $\mu\text{m}^2$ )	12811307	4110861
Total Stain Area ( $\mu\text{m}^2$ )	95351.09	62105
Average Stain Intensity	158.667	191
Mean Vessel Area ( $\mu\text{m}^2$ )	143.066	213
Median Vessel Area ( $\mu\text{m}^2$ )	99.6743	121
Standard Deviation of Vessel Area	135.16	302
Mean Vessel Perimeter ( $\mu\text{m}$ )	65.1583	78
Median Vessel Perimeter ( $\mu\text{m}$ )	50.17	58
Standard Deviation of Vessel Perimeter ( $\mu\text{m}$ )	44.7857	66
Mean Lumen Area ( $\mu\text{m}^2$ )	16.9771	46
Median Lumen Area ( $\mu\text{m}^2$ )	0	0
Standard Deviation of Lumen Area ( $\mu\text{m}^2$ )	65.3677	196
Mean Vascular Area ( $\mu\text{m}^2$ )	114.144	157
Median Vascular Area ( $\mu\text{m}^2$ )	82.3068	98
Standard Deviation of Vascular Area ( $\mu\text{m}^2$ )	105.778	195
Mean Vessel Wall Thickness ( $\mu\text{m}$ )	1.5273	1.66
Median Vessel Wall Thickness ( $\mu\text{m}$ )	1.47374	1.55
Standard Deviation of Vessel Wall Thickness ( $\mu\text{m}$ )	0.402804	0.53




## Saving and Using the Algorithm Macro

To use the algorithm for analysis, you will need to save the algorithm macro, which consists of all of the settings you have adjusted.

If you are satisfied with the results based on your parameter changes, save the macro. You are asked to enter a name for the macro; type a name that will help you identify the macro in the future. Because the algorithm has been optimized for a particular specimen, we suggest you use a macro name that will remind you which specimen was analyzed. Now, when you bring up the Analysis window, you will be able to select your macro to run an analysis. See the *Aperio Image Analysis User's Guide* for information.

## Example - Oncology Research with Xenografts

The following example shows using the microvessel analysis algorithm to quantify vascular changes extending outward from a xenograft, and then comparing with normal tissue vasculature.



## For More Information

For details on using the Aperio products mentioned in this guide, see:

- ▶ *eSlide Manager Operator's Guide* and *eSlide Manager Administrator's Guide* for information on using eSlide Manager.
- ▶ *ImageScope User's Guide* for details on using the ImageScope viewer.
- ▶ User's Guide for your Aperio scanner for information on scanning glass slides.

# Index

## A

analysis 25  
modes 26

## B

background stain 12

## C

color deconvolution 16  
completing vessels  
  measuring gaps 22  
  Region Joining Parameter 21  
  Vessel Completion Parameter 22

## D

dark staining threshold 19  
deconvolution 18  
definitions 14

## E

endothelial cell 14  
endothelial stain 11

## H

histograms 11  
  bins 11  
  creating 26  
  values 11

## I

ImageScope drawing tools 17  
input parameters 10  
installation 8

intended use 8

## L

light staining threshold 20  
lumen 14

## M

mark-up image type 11, 17  
  tuning 17  
metrics 25  
microvessel analysis algorithm 7  
monitor requirements 8

## P

parameters 10

## R

region 14  
Region Joining Parameter 11, 21  
results 12

## S

selecting the algorithm 15  
stain color vectors 17  
staining thresholds 11, 19














## V

vascular area 14  
vessel 14  
Vessel Completion Parameter 11, 22



## Symbols

- The following symbols may appear on your product label or in this user's guide:

	Manufacturer
	Date of manufacture (year - month - day)
	European Union Authorized Representative
	In vitro diagnostic device
	Serial number
	Catalog number
	Relative humidity range
	Storage temperature range
	Electronic and electrical equipment waste disposal
	The exclamation point within an equilateral triangle is intended to alert you to the presence of important operating and maintenance (servicing) instructions. <i>Le point d'exclamation dans un triangle équilatéral vise à avertir l'utilisateur qu'il s'agit d'instructions d'utilisation et d'entretien importantes.</i>
	The lightning flash with arrowhead symbol within an equilateral triangle is intended to alert you to the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons. <i>Le symbole de l'éclair avec la pointe de flèche dans un triangle équilatéral vise à avertir l'utilisateur que le boîtier du produit présente une « tension dangereuse » non isolée d'une amplitude suffisante pour constituer un risque d'électrocution.</i>
	The flat surface with waves symbol within an equilateral triangle is intended to alert you to the presence of hot surfaces which could cause burn damage. <i>Le symbole d'une surface plane et de vagues dans un triangle équilatéral vise à avertir l'utilisateur de la présence de surfaces chaudes qui peuvent causer des brûlures.</i>
	The UV lamp within an equilateral triangle is intended to alert you to the presence of UV light within the product's enclosure that may be of sufficient magnitude to constitute a risk to the operator. <i>La lampe UV dans un triangle équilatéral vise à avertir l'utilisateur de la présence de rayonnement UV dans le boîtier du produit qui peut être d'une amplitude suffisante pour constituer un risque pour l'utilisateur.</i>

