

BLOOD GLUCOSE UPLOADING FROM CGM'S AND METERS

IN ORDER TO PROVIDE THE BEST CARE FOR YOUR CHILD, PLEASE UPLOAD ALL DIABETES DEVICES 3 DAYS PRIOR TO YOUR CHILD'S APPOINTMENT OR IF YOU ARE CALLING, MYCHART OR EMAILING FOR A BG REVIEW PLEASE UPLOAD INFORMATION PRIOR.

Glooko (Meters, Omnipod Insulin Pump)

<https://join.glooko.com>

- Please check if your meter and/or insulin pump is compatible with Glooko at <https://www.glooko.com/compatibility>
- If your device is compatible, please go to Glooko's website and either create an account or click "Log In" if you already have a Glooko account. After logging into your account, click "settings>account". Click "add new code" and then "continue". Our ProConnect code is **wupeds**, then click "submit".
- Please call/email our office or send us a message through MyChart after uploading.

Tandem T:SLIM X2 Pump with CONTROL-IQ Technology

<https://tconnect.tandemdiabetes.com>

- Log into t:connect.
- Upload your t:slim X2 pump.
- Please call/email our office or send a message through MyChart after uploading.
- Please provide us with your t:connect log in information (email & password) so that we can review your upload.

Medtronic CareLink (Medtronic 630G/670G/770G insulin pump only)

<http://medtronicdiabetes.com/carelink>

- Please either click "Sign Up Now" or log in if you already have an account.
- Please call/email our office or send us a message through MyChart after uploading.

Clarity (Dexcom G5/G6 CGM only)

<https://clarity.dexcom.com>

- If you use your iPhone/Android as your Dexcom receiver, please call our office and provide an email address we can send an invitation to for sharing.
- If you use a Dexcom G5/G6 receiver, please upload your receiver to the Clarity website and provide us with an email address we can send an invitation to for sharing.
- Please call/email our office or send us a message through MyChart after uploading.

LibreView (FreeStyle Libre CGM only)

<https://www.libreview.com>

- Please log in or create an account and upload your device.
- Please call/email our office or send us a message through MyChart after uploading and provide us your email address so that we can email you an invitation for sharing data.

One Touch Verio Meters

<https://onetouchreveal.com/login>

- Please load Reveal APP on your mobile phone
- Please log in or create an account
- Clinic Code: **JD QE TQ**
- Select the meter; Reflect, Flex and/or Sync
- Follow App instruction to pair with meter

Unable To Upload Your Device

- Please MyChart, email, or fax us two weeks of blood sugars to our office.
 - Can also send us pictures of your blood sugar log.
- **Please include patient's name and date of birth.**

Insulin Pump Representative Information

Tandem(T-Slim):

Sales - Melody Fox

mfox@tandemdiabetes.com

Cell: 314.449.4771

Specialist/Trainer – Sarah Seaman

sseaman@tandemdiabetes.com

Cell: 314.200.4664

IT Support Number: 877.801.6901

Medtronic:

Sales – Nick Harper

nicholas.s.harper@medtronic.com

Cell 314.915.4444

Specialist/Trainer – Sara Pryor

sara.pryor@medtronic.com

Cell: 314.625.6905

IT Support Number: 800.284.4416

Omnipod:

Sales – Michaela Esch

mesch@insulet.com

Cell: 314.412.4323

Specialist/Trainer – Meah Simmons

Cell: 618.402.6769

Please call/email/MyChart us with any questions!

Washington University

Pediatric Endocrinology & Diabetes

Phone: (314) 454-6051

Fax: (314) 454-6225

Email: pedsdiabetesteam@email.wustl.edu

MyChart: <https://www.mypatientchart.org>